



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet rack for 16 cryoboxes 136x136x53 mm

€145.10
plus VAT &
Shipping

Product Images



Description

The **Cryomaster® cabinet racks** for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. A practical folding handle is located on the front and back of the cabinet rack for easy removal.

Additional Information

| | |
|----------------------------------|---|
| No. | 4-6203 |
| Manufacturer (Brand) | Cryomaster |
| EAN | 4058072160685 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Wide | 139 mm |
| Height | 224 mm |
| Length | 569 mm |
| Depth | 569 mm |
| Number of compartments in height | 4 |
| Number of compartments in depth | 4 |
| suitable for cooling in | Ultra Deep Freezer Freezer Refrigerator |
| Type slide-in | Cabinet racks for 50 boxes |

