



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet rack for 20 cryoboxes 136x136x53 mm

€168.80
plus VAT &
Shipping

Product Images



Description

The **Cryomaster® cabinet racks** for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. A practical folding handle is located on the front and back of the cabinet rack for easy removal.

Additional Information

No.	4-6204
Manufacturer (Brand)	Cryomaster
EAN	4058072158088
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	139 mm
Height	279 mm
Length	569 mm
Depth	565 mm
Number of compartments in height	5
Number of compartments in depth	4
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator
Type slide-in	Cabinet racks for 50 boxes

