



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® cabinet rack narrow for 4 cryoboxes 136x136x275mm

**€168.80**  
**plus VAT &**  
**Shipping**

### Product Images



Abbildung ähnlich

## Description

---

The **Cryomaster® cabinet racks** for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. A practical folding handle is located on the front and back of the cabinet rack for easy removal.

## Additional Information

---

No.	4-6207
Manufacturer (Brand)	Cryomaster
EAN	4058072157753
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	139 mm
Height	279 mm
Length	569 mm
Depth	569 mm
Number of compartments in height	1
Number of fans length	4
Capacity	4
Type slide-in	Cabinet racks for 50 boxes
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator

