




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® cabinet box outer dimensions WxDxH=278x685x279 mm

**€233.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Cryomaster **cabinet box** is an alternative to the usual cabinet drawers or racks.

It is used to store samples or storage containers of different sizes.

Suitable for use in refrigerators or freezers.

A wide grip edge on the front allows easy removal of the box.

The box is closed on 3 sides; access is possible through the open top.

## Additional Information

---

No.	4-6240
Manufacturer (Brand)	Cryomaster
EAN	4058072191658
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	278 mm
Height	279 mm
Length	685 mm
Depth	685 mm
Type slide-in	Cabinet boxes
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator

