



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet box outer dimensions WxDxH=278x705x224 mm

€255.80
plus VAT &
Shipping

Product Images

Description

The Cryomaster **cabinet box** is an alternative to the usual cabinet drawers or racks.

It is used to store samples or storage containers of different sizes.

Suitable for use in refrigerators or freezers.

A wide grip edge on the front allows easy removal of the box.

The box is closed on 3 sides; access is possible through the open top.

Additional Information

No.	4-6254
Manufacturer (Brand)	Cryomaster
EAN	4058072191795
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	278 mm
Height	224 mm
Length	705 mm
Depth	705 mm
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator
Type slide-in	Cabinet boxes

