



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® cabinet rack incl. 16 cardboard cryoboxes without grid 133x133x50mm, 4 deep, 4 high

Special Price

€62.00 was

€128.10

plus VAT &  
Shipping

### Product Images



## Description

---

The Cryomaster® cabinet racks for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. A practical folding handle is located on the front and back of each cabinet rack for easy removal. The compartment dimensions are suitable for the internationally common box dimensions of 133 x 133 x 50 mm (L x W x H). The capacity of this cabinet rack is 16 boxes (4 deep, 4 high).

## Additional Information

No.	4-6007
Manufacturer (Brand)	Cryomaster
EAN	4058072087081
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	544 mm
Width/inside	138 mm
Height	220 mm
Height inside	55.0 mm
Length	140 mm
Depth	140 mm
Length inside	136 mm
Number of compartments in height	4
Number of compartments in depth	4
Number of fans length	1
suitable for cooling in	Ultra Deep Freezer / Deep Freezer
Maximum width box	133 mm
Maximum height box	50 mm
Maximum length box	133 mm
Type slide-in	Cabinet racks for 50 boxes

