



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet rack f. 50erBox, 3x4 compartments, 142 x 426 x 224 mm

Special Price

~~€99.00~~ was

€148.50

**plus VAT &
Shipping**

Product Images



Description

For all boxes with external dimensions 136 x 136 x 53 mm (W x D x H).

Suitable e.g. for neolab boxes 2-2695 ff, 2-1916, 2-2900 or 7-8030 ff.

Basic dimension: 142 mm width

Additional Information

No.	2-2742
Manufacturer (Brand)	Cryomaster
EAN	4058072068790
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	426 mm
Width/inside	140 mm
Height	224 mm
Height inside	55.0 mm
Length	142 mm
Depth	142 mm
suitable for	Boxes/boxes
Length inside	140 mm
Racks properties	stackable / interlockable
Number of compartments in height	4
Number of compartments in depth	3
Number of fans length	1
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container
Maximum width box	136 mm
Maximum height box	53 mm
Maximum length box	136 mm
Type slide-in	Cabinet racks for 50 boxes

