



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet rack f. neoBox-50, 4x5 compartments, 146 x 362 x 299 mm

€119.00
plus VAT &
Shipping

Product Images



Description

For all cryo boxes with external dimensions 140 x 86 x 55 mm (W x D x H).

Specially manufactured for the neoBox-50, order no. 2-1966 and 2-1941 ff.

Basic dimension: 146 mm width

Additional Information

No.	2-8269
Manufacturer (Brand)	Cryomaster
EAN	4058072073985
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	146 mm
Width/inside	144 mm
Height	299 mm
Height inside	60.0 mm
Length	362 mm
Depth	362 mm
Length inside	88 mm
Number of compartments in height	5
Number of compartments in depth	1
Number of fans length	4
Maximum width box	88 mm
Maximum height box	55 mm
Maximum length box	140 mm
Type slide-in	Cabinet racks for neoBox-50
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container

