



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® chest rack incl. 11 cardboard cryoboxes without grid 133x133x50mm

Special Price

€52.00 was

€92.00

plus VAT &  
Shipping

### Product Images



## Description

---

The Cryomaster® chest racks for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. The integrated retaining rod prevents the boxes from falling out of the chest rack. A convenient folding handle is located on the top of the chest rack for easy removal. The compartment dimensions are suitable for the internationally common box dimensions of 133 x 133 x 50 mm (L x W x H).

The capacity of this chest rack is 11 boxes.

## Additional Information

No.	4-6003
Manufacturer (Brand)	Cryomaster
EAN	4058072087043
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	silver
Material	Stainless steel
sterile	No
Wide	140 mm
Width/inside	138 mm
Height	605 mm
Height inside	55.0 mm
Length	140 mm
Depth	140 mm
Length inside	138 mm
Number of compartments in height	11
Number of compartments in depth	1
Number of fans length	1
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Maximum width box	135 mm
Maximum height box	50 mm
Maximum length box	135 mm
Type slide-in	Chest racks for 50 boxes

