



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® chest rack incl. 12 cardboard cryoboxes without grid 133x133x50mm

Special Price

€55.00 was

€94.00

**plus VAT &
Shipping**

Product Images

Description

The Cryomaster® chest racks for sample storage in freezers.

Precise manufacturing from high-quality stainless steel enables optimum use of space with maximum stability in your freezer. The integrated retaining rod prevents the boxes from falling out of the chest rack. A convenient folding handle is located on the top of the chest rack for easy removal. The compartment dimensions are suitable for the internationally common box dimensions of 133 x 133 x 50 mm (L x W x H).

The capacity of this chest rack is 12 boxes.

Additional Information

No.	4-6004
Manufacturer (Brand)	Cryomaster
EAN	4058072087050
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	140 mm
Width/inside	138 mm
Height	660 mm
Height inside	55.0 mm
Length	140 mm
Depth	140 mm
Length inside	138 mm
Number of compartments in height	12
Number of compartments in depth	1
Number of fans length	1
Maximum width box	135 mm
Maximum height box	50 mm
Maximum length box	135 mm
Type slide-in	Chest racks for 50 boxes
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

