



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® chest rack with safety bar for 11 cryoboxes 136x136x53 mm

Special Price

€65.00 was

€101.60

**plus VAT &
Shipping**



Description

Cryomaster® stainless steel chest racks for all boxes with external dimensions up to 136 x 136 x 53 mm (W x D x H). Specially manufactured for the neoLab cryo box (order no. 2-2695 ff.). Also suitable for order no. 2-1916, 2-2700, 2-2900 or 7-8030 ff.

Additional Information

No.	4-6205
Manufacturer (Brand)	Cryomaster
EAN	4058072156923
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	141 mm
Height	611 mm
Length	140 mm
Depth	148 mm
Number of compartments in height	11
Maximum width box	135 mm
Maximum height box	50 mm
Maximum length box	135 mm
Type slide-in	Chest racks for 50 boxes
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator

