



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® chest rack f. 100 box, 6 compartments, 141 x 141 x 636 mm

€97.00
plus VAT &
Shipping



Description

For all cryo boxes with external dimensions up to 136 x 136 x 103 mm (W x D x H).

Specially manufactured for the neoLab cryo box (order no. 2-2702).
Also suitable for order no. 2-2902 or 7-8020 ff.

Basic dimensions: 141 x 141 mm (W x D)

Do you require a different number of compartments than specified? No problem, just ask us!

Additional Information

No.	2-2715
Manufacturer (Brand)	Cryomaster
EAN	4058072068677
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	141 mm
Width/inside	140 mm
Height	636 mm
Height inside	105.0 mm
Length	141 mm
Depth	141 mm
suitable for	Boxes/boxes
Length inside	140 mm
Number of compartments in height	6
Number of compartments in depth	1
Number of fans length	1
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container
Maximum width box	136 mm
Maximum height box	103 mm
Maximum length box	136 mm
Type slide-in	Chest racks for 100 boxes

