



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Cryomaster® chest rack f. 50 box, 8
compartments, 141 x 141 x 446 mm**

**€78.00
plus VAT &
Shipping**



Description

For all boxes with external dimensions up to 136 x 136 x 53 mm (W x D x H).

Specially manufactured for the neolab cryogenic cardboard box (order no. 2-2695 ff.).
Also suitable for order no. 2-1916, 2-2700, 2-2900 or 7-8030 ff.

Basic dimensions: 141 x 141 mm (W x D)

Additional Information

| | |
|----------------------------------|---|
| No. | 2-8513 |
| Manufacturer (Brand) | Cryomaster |
| EAN | 4058072010409 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Wide | 141 mm |
| Width/inside | 140 mm |
| Height | 446 mm |
| Height inside | 58.0 mm |
| Length | 141 mm |
| Depth | 141 mm |
| suitable for | Boxes/boxes |
| Length inside | 140 mm |
| Number of compartments in height | 8 |
| Number of compartments in depth | 1 |
| Number of fans length | 1 |
| suitable for cooling in | Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container |
| Maximum width box | 136 mm |
| Maximum height box | 53 mm |
| Maximum length box | 136 mm |
| Type slide-in | Chest racks for 50 boxes |

