



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® chest rack f. neoBox-81, 5 compartments, 138 x 138 x 281 mm

€93.60
plus VAT &
Shipping



Description

For all boxes with external dimensions up to 133 x 133 x 53 mm (W x D x H).

Specially manufactured for the neoBox-81 (order no. 2-2910 ff).

Also suitable for order no. 2-2580 ff or 7-8041 ff.

Basic dimensions: 138 x 138 mm (W x D)

Additional Information

No.	2-2716
Manufacturer (Brand)	Cryomaster
EAN	4058072068684
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	138 mm
Width/inside	137 mm
Height	281 mm
Height inside	56.0 mm
Length	138 mm
Depth	138 mm
suitable for	Boxes/boxes
Number of compartments in height	5
Number of compartments in depth	1
Number of fans length	1
Capacity	5
Length inside	137 mm
Maximum width box	133 mm
Maximum height box	53 mm
Maximum length box	133 mm
suitable for cooling in	Nitrogen tank

