



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoStackRack® Rack for stacking, white

**€6.70** As low as

**€5.89**

**plus VAT &  
Shipping**

### Product Images



## Description

---

### **Individual element for** neoStapelRack®

Suitable for common standard cryogenic storage boxes with heights of 50-53 mm.

There are large openings on the side walls and bottom for better cold or heat transfer.

Each stacking platform has a crescent-shaped recess at the front for easier removal of the storage boxes.

### **Dimensions basic element outside:**

146 x 150 x 58 mm (W x D x H)

### **Dimensions basic element inside:**

140 x 147 x 55 mm

**Material:** PC white, autoclavable

Temperature resistant from -130°C to +130°C

## Additional Information

---

No.	2-2999
Manufacturer (Brand)	neoStapelRack
EAN	4058072000110
Transport temperature	Room temperature
Color	White
Material	Polycarbonate (PC)
Wide	150 mm
Width/inside	140 mm
Height	58 mm
Height inside	55.0 mm
Length	146 mm
Depth	146 mm
suitable for	Boxes/boxes
Length inside	147 mm
Racks properties	stackable / interlockable
TBST MAX	130 °C
TBST MIN	-130 °C
Number of compartments in height	1
Number of compartments in depth	1
Number of fans length	1
suitable for cooling in	Ultra Deep Freezer / Deep Freezer
Maximum width box	140 mm
Maximum height box	53 mm
Maximum length box	147 mm

