




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 4-sided pinball rack size 1, blue

**€27.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

All sides of the rack can be used individually.

autoclavable, temperature resistant from -90°C to +121°C.

Material: Polypropylene

Dimensions: 197 x 95 x 115 mm (L x D x H)

### Flipper rack 1

8 holes 30 mm Ø

20 holes 20 mm Ø

24 holes 17 mm Ø

40 holes 13 mm Ø

## Additional Information

No.	2-2773
Manufacturer (Brand)	neoLab
EAN	4058072069001
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Wide	95 mm
Height	115 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
for vessels with volume up to	50 ml 30 ml 15 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	13.0 mm 17.0 mm 20.0 mm 30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels centrifuge tubes test tubes/tubes

