



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Adapt-a-rack™, POM, blue

**€18.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

neoLab-Adapt-a-Rack™ - a rack for different tube sizes. The flexible rack is suitable for 4 tubes with a diameter of 12 to 30 mm (5 ml - 50 ml). The special holder keeps the tubes absolutely straight.

The individual racks can be connected to each other. Due to the open design the labeling of the tubes can be read very well.

Rack size: 18 x 5.5 x 7.5 cm (LxWxH)

## Additional Information

---

No.	2-6627
Manufacturer (Brand)	neoLab
EAN	4058072008239
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polyoxymethylene (POM)
Wide	55 mm
Height	75 mm
Racks features	for different vessel types/volumes
Number of floors	1
Capacity	4
Squares (WxL)	2 x 4
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	12.0 mm 30.0 mm
for vessels with volume up to	5 ml 50 ml
for micro vessel/tube shape	for all shapes
suitable for cooling in	Refrigerator

