



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Adapt-a-rack™, POM, yellow

€18.00
plus VAT &
Shipping

Product Images



Description

neoLab-Adapt-a-Rack™ - a rack for different tube sizes. The flexible rack is suitable for 4 tubes with a diameter of 12 to 30 mm (5 ml - 50 ml). The special holder keeps the tubes absolutely straight.

The individual racks can be connected to each other. Due to the open design the labeling of the tubes can be read very well.

Rack size: 18 x 5.5 x 7.5 cm (LxWxH)

Additional Information

No.	2-6630
Manufacturer (Brand)	neoLab
EAN	4058072008260
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polyoxymethylene (POM)
Wide	55 mm
Height	75 mm
Racks features	for different vessel types/volumes
Number of floors	1
for vessels with volume up to	5 ml 50 ml
suitable for cooling in	Refrigerator
Squares (WxL)	2 x 4
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	12.0 mm 30.0 mm
for micro vessel/tube shape	for all shapes

