



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Combi Tube Rack, yellow

€18.40
plus VAT &
Shipping

Product Images



Description

All 4 sides of the rack can be used for different tubes.

Each compact unit includes:

4 holes 30 mm Ø for 50 ml tubes

12 holes 17 mm Ø for 15 ml centrifuge tubes

32 holes 12 mm Ø for 1.5 and 2 ml reaction tubes

32 holes 8 mm Ø for 0.5 ml reaction tubes

Dimensions: 180 x 90 x 52 mm (L x H x D)

Material: Polypropylene

Temperature resistant from -90°C to +121°C

Autoclavable

Additional Information

No.	2-2686
Manufacturer (Brand)	neoLabLine
EAN	4058072068455
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	52 mm
Height	90 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
Capacity	variable (depending on vessel volume)
Squares (WxL)	variable
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	10.0 mm 12.0 mm 20.0 mm 30.0 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml 15 ml 50 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels Centrifuge tubes
suitable for cooling in	Freezer/ freezer

