



neoLabLine Combi Tube Rack, blue

Product Images



Description

All 4 sides of the rack can be used for different tubes.

Each compact unit includes:

4 holes 30 mm Ø for 50 ml tubes

12 holes 17 mm Ø for 15 ml centrifuge tubes

32 holes 12 mm Ø for 1.5 and 2 ml reaction tubes

32 holes 8 mm Ø for 0.5 ml reaction tubes

Dimensions: 180 x 90 x 52 mm (L x H x D)

Material: Polypropylene

Temperature resistant from -90°C to +121°C

Autoclavable

Additional Information

No.	2-2688
Manufacturer (Brand)	neoLabLine
EAN	4058072068479
Color	blau
Material	Polypropylen (PP)
Wide	52 mm
Height	90 mm
suitable for	variable Anwendungen
Racks features	stapel-/zusammensteckbar mit besonderer chemikalischer Beständigkeit inkubatorgeeignet für verschiedene Gefäßtypen/Volumina alphanumerisch codiert
Number of floors	1
for vessels with volume up to	0.5 ml 1.5 ml 2 ml 15 ml 50 ml
suitable for cooling in	Gefrierschrank/ -truhe
Squares (WxL)	variabel
Type frame	Kombi-Rack für verschiedene Gefäße und Volumina
for vessels with DM up to	10.0 mm 12.0 mm 20.0 mm 30.0 mm
for micro vessel/tube shape	für alle Formen
for type (micro) vessel	Reaktionsgefäße Zentrifugenröhrchen

