



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Combi Tube Rack Extra Pack, all colors, 6 pcs/pack

€72.50
plus VAT &
Shipping

Product Images



Description

All 4 sides of the rack can be used for different tubes.

Each compact unit includes:

4 holes 30 mm Ø for 50 ml tubes

12 holes 17 mm Ø for 15 ml centrifuge tubes

32 holes 12 mm Ø for 1.5 and 2 ml reaction tubes

32 holes 8 mm Ø for 0.5 ml reaction tubes

Dimensions: 180 x 90 x 52 mm (L x H x D)

Material: Polypropylene

Temperature resistant from -90°C to +121°C

Autoclavable

Order No. 2-2693

1x each: transparent, yellow, green, blue, pink, orange

Additional Information

No.	PK-22693
Manufacturer (Brand)	neoLabLine
EAN	4058072068509
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky) yellow green blue pink orange
Material	Polypropylene (PP)
Wide	52 mm
Height	90 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
for vessels with volume up to	0.5 ml 1.5 ml 2 ml 15 ml 50 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	10.0 mm 12.0 mm 20.0 mm 30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels Centrifuge tubes

