




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Wire rack for tubes 16 mm, 108 places, red

**€37.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Autoclavable red wire rack with epoxy coating for tubes up to 16 mm in diameter.

9x 12 square holes for test tubes, centrifuge tubes or other vessels.

Reinforced corners, additional grid bars at the bottom and a second grid at mid-height provide stability and a secure hold.

## Additional Information

---

No.	2-5089
Manufacturer (Brand)	neoLab
EAN	4058072004286
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Steel epoxy coated
Wide	188 mm
Height	65 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	16 x 16 mm
Squares (WxL)	9 x 12
Type frame	Test tube rack
for vessels with DM up to	16.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes

