



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand test tube racks, PP, for 84 tubes up to diameter 13 mm, red

**€17.49** As low as  
**€13.45**  
plus VAT &  
Shipping

### Product Images



## Description

---

Red frame made of polypropylene, stable even in moving water.  
Alphanumeric row coding for better identification of the glasses.  
The racks are supplied folded and are easy to put together.  
Stackable

## Additional Information

No.	BN-1255
Manufacturer (Brand)	Brand
old neoLab article no.	6-2510
EAN	4058072206345
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Wide	126 mm
Height	75 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
suitable for cooling in	Freezer/ freezer
Tray size	13 x 13 mm
Squares (WxL)	6 x 14
Type frame	Test tube rack
for vessels with DM up to	13.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes

