



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand test tube racks, PP, for 40 tubes up to diameter 20 mm, blue

**€17.49** As low as  
**€13.45**  
plus VAT &  
Shipping

### Product Images



## Description

---

Blue polypropylene frame, stable even in moving water.  
Alphanumeric row coding for better identification of the glasses.  
The racks are delivered folded and are easy to put together.  
Stackable

## Additional Information

No.	BN-1262
Manufacturer (Brand)	Brand
old neoLab article no.	6-2522
EAN	4058072206390
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Wide	126 mm
Height	75 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
Tray size	20 x 20 mm
Squares (WxL)	4 x 10
Type frame	Test tube rack
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes
suitable for cooling in	Freezer/ freezer

