



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand test tube racks, PP, for 21 tubes up to diameter 30 mm, blue

€17.49 As low as
€13.45
plus VAT &
Shipping

Product Images



Description

Blue polypropylene frame, stable even in moving water.
Alphanumeric row coding for better identification of the glasses.
The racks are delivered folded and are easy to put together.
Stackable

Additional Information

No.	BN-1269
Manufacturer (Brand)	Brand
old neoLab article no.	6-2524
EAN	4058072206451
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Wide	126 mm
Height	88 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
Tray size	30 x 30 mm
Squares (WxL)	3 x 7
Type frame	Test tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes
suitable for cooling in	Freezer/ freezer

