



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand test tube racks, PP, for 55 tubes up to diameter 16 mm, red

€17.49 As low as
€13.45
plus VAT &
Shipping

Product Images



Description

Red frame made of polypropylene, stable even in moving water.
Alphanumeric row coding for better identification of the glasses.
The racks are delivered folded and are easy to put together.
Stackable

Additional Information

No.	BN-1273
Manufacturer (Brand)	Brand
old neoLab article no.	6-2511
EAN	4058072206499
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Wide	126 mm
Height	75 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
suitable for cooling in	Freezer/ freezer
Tray size	16 x 16 mm
Squares (WxL)	5 x 11
Type frame	Test tube rack
for vessels with DM up to	16.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes

