



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Scintillation Vessel Rack, 24 places 30 mm Ø

€62.90
plus VAT &
Shipping

Product Images



Description

Stable rack with 24 openings for easy and convenient handling of scintillation vessels, for storage, transport or sample preparation. Made of chemical-resistant polypropylene, autoclavable up to 121°C.
Stable design with 2 carrying handles.

Intermediate shelf height: 44 mm

Holes:

3 rows with 8 holes of 30 mm Ø each

Additional Information

No.	1-3107
Manufacturer (Brand)	neoLab
EAN	4058072035020
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	125 mm
Height	93 mm
suitable for	Vials/ Vials
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
DM Bore	30 mm
for vessels with volume up to	20 ml
suitable for cooling in	Freezer/ freezer
Tray size	30 x 30 mm
Squares (WxL)	3 x 8
Type frame	Scintillation Vessel Rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	round conical. with stand ring
for type vial/vial	Vial. 30 mm Ø

