



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Combi-Cube Rack for Tubes, 5 pcs/pack

€67.60
plus VAT &
Shipping

Product Images



Description

One for all!

This versatile rack made of PP can be reassembled as needed.

It accommodates 4 x 50 ml tubes, 10 x 15 ml conical tubes, 12 x 12 mm tubes (75 or 100 mm long) or 16 x 1.5 ml to 2.0 ml reaction tubes.

Autoclavable up to 121°C.

Pack of 5 assorted color racks in blue, green, yellow, orange and pink.

Dimensions: 11 x 11 x 11 cm

Additional Information

No.	2-1829
Manufacturer (Brand)	neoLab
EAN	4058072058111
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue yellow green orange pink
Material	Polypropylene (PP)
Wide	110 mm
Height	110 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
Squares (WxL)	variable
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	11.0 mm 12.0 mm 16.0 mm 30.0 mm
for vessels with volume up to	1.5 ml 2 ml 15 ml 50 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels Centrifuge tubes
suitable for cooling in	Freezer/ freezer

