



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® LaboBox system: rack for vials up to
12.5 mm, 5 x 10 places, transp.**

€9.76 As low as

€8.92

**plus VAT &
Shipping**

Product Images



Description

The core of the LaboBox® system are the racks for reaction vessels, vials, cryotubes or vials, which you can also use as table racks. The base is closed.

For a better overview, the compartments are alphanumerically labeled.

The rack has two carrying straps on the side.

The racks are temperature-resistant from -130°C to +130°C and autoclavable.

Material: PC

Dimensions: 164 x 76 x 28 mm (L x W x H)

Order No. 2-2312 Rack for vials up to 12.5 mm Ø

e.g. cryo vials 1.2 ml (3-6000), micro vials 1.5 ml (2-1000), sample vials 2.0 ml (2-7220 /2-7230)

Additional Information

No.	2-2312
Manufacturer (Brand)	neoLab
EAN	4058072060510
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polycarbonate (PC)
Wide	76 mm
Height	28 mm
suitable for	Vials/ Vials
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	12.5 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Squares (WxL)	5 x 10
Type frame	Vial/Vial Rack
for vessels with DM up to	12.5 mm
for micro vessel/tube shape	conical. with stand ring
for type (micro) vessel	Cryotube
for type vial/vial	Vial. 12.5 mm Ø

