



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® LaboBox system: rack for vials up to 15 mm, 4 x 9 places, transp.

€9.76 As low as

€8.92

**plus VAT &
Shipping**

Product Images



Description

The core of the LaboBox® system are the racks for reaction vessels, vials, cryotubes or vials, which you can also use as table racks. The base is closed.

For a better overview, the compartments are alphanumerically labeled.

The rack has two carrying straps on the side.

The racks are temperature-resistant from -130°C to +130°C and autoclavable.

Material: PC

Dimensions: 164 x 76 x 28 mm (L x W x H)

Order No. 2-2313 Rack for vials up to 15 mm Ø

e.g. micro vials 4.0 ml (2-1001), sample vials 4.0 ml (2-7221 / 2-7231)

Additional Information

| | |
|-------------------------------|--|
| No. | 2-2313 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072060527 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Polycarbonate (PC) |
| Wide | 76 mm |
| Height | 28 mm |
| suitable for | Vials/ Vials |
| Racks features | with special chemical resistance suitable for incubators for different vessel types/volumes alphanumerically coded |
| Number of floors | 1 |
| DM Bore | 15 mm |
| for vessels with volume up to | 4 ml |
| suitable for cooling in | Ultra Deep Freezer/ Deep Freezer |
| Squares (WxL) | 4 x 9 |
| Type frame | Vial/Vial Rack |
| for vessels with DM up to | 15.0 mm |
| for micro vessel/tube shape | conical. with stand ring |
| for type vial/vial | Vial. 15 mm Ø |

