



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine multifunctional rack in basket form, yellow

€12.50
plus VAT &
Shipping

Product Images



Description

neoLabLine **multifunctional rack in basket form**

Practical basket as transport or storage option in the laboratory.
Consists of a base plate, 2 holding frames and a handle.

- Stackable
- Autoclavable

Additional Information

No.	6-2611
Manufacturer (Brand)	neoLabLine
EAN	4058072156022
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	120 mm
Height	70 mm
Depth	230 mm
Racks features	for different vessel types/volumes alphanumerically coded
for vessels with volume up to	50 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	2 x 4
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes

