



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Multismart Rack, white

**€18.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

Just about anything fits into the **neolabLine Multismart Rack** made of chemically resistant polypropylene. Simply insert the pins into the plate according to individual requirements and thus store or transport various containers and vessels.

Test tubes, plates, vials or other vessels with different diameters find a firm and stable hold in the rack. The Multismart Rack can also be used as a draining rack

- With 105 pins for individual assembly
- Autoclavable

## Additional Information

No.	2-1199
Manufacturer (Brand)	neolabLine
EAN	4058072155759
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	100 mm
Height	65 mm
Depth	230 mm
Racks features	for different vessel types/volumes
for vessels with volume up to	10 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	7 x 17
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	12.0 mm
for micro vessel/tube shape	for all shapes

