



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Test tube racks for autosampler, white, 4 x 10 hole for tubes with 20 mm Ø

**€36.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Stackable test tube racks for autosamplers in various designs

The easily stackable racks are made of robust polypropylene and can therefore be autoclaved without any problems. They fit into all common autosamplers (corresponding size) and are equally stable even in moving water. Alphanumeric row coding for better identification of the tubes.

## Additional Information

---

No.	6-2542
Manufacturer (Brand)	neoLab
EAN	4058072156640
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Wide	104 mm
Height	64 mm
Racks features	incubator suitable for different vessel types/volumes alphanumerically coded
DM Bore	20.5 mm
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Squares (WxL)	4 x 10
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes

