



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Test tube rack for autosampler, red, 3 x 8 holes for tubes with 25 mm Ø

**€36.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Stackable test tube racks for autosamplers in various designs

The easily stackable racks are made of robust polypropylene and can therefore be autoclaved without any problems. They fit into all common autosamplers (corresponding size) and are equally stable even in moving water. Alphanumeric row coding for better identification of the tubes.

## Additional Information

---

No.	6-2548
Manufacturer (Brand)	neoLab
EAN	4058072161262
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Wide	104 mm
Height	64 mm
Racks features	stackable/collapsible suitable for incubator for different vessel types/volumes alphanumerically coded
DM Bore	25.5 mm
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Squares (WxL)	3 x 8
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	25.0 mm
for micro vessel/tube shape	for all shapes

