



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rota-Rack for reaction tubes

€61.90
plus VAT &
Shipping

Product Images



Description

Colorful and extremely practical, the Rota-Rack® comes in two versions for different tube sizes. The four colored elements (blue, green, pink, yellow) are connected to each other, but each element can be rotated in the desired direction so that tubes of different sizes can be stored and transported in one rack. Autoclavable.

Rota-Rack® for reaction tubes

6 x 15 ml

9 x 1.5/2.0 ml

12 x 0.5/0.6 ml

32 x 0.2 PCR tubes / 4 x 8-strips

Dimensions: 208 x 100 x 100 mm

Additional Information

No.	2-6616
Manufacturer (Brand)	neoLab
EAN	4058072008123
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue yellow green pink
Material	Polypropylene (PP)
Wide	100 mm
Height	100 mm
suitable for	variable applications
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
Squares (WxL)	variable
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes Combi rack for different tubes and volumes
for vessels with volume up to	0.2 ml 0.5 ml 0.6 ml 1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	PCR Reaction Vessels Reaction Vessels Centrifuge Tubes
suitable for cooling in	Freezer/ freezer

