



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rota-Rack for large tubes

€63.70
plus VAT &
Shipping

Product Images



Description

Colorful and extremely practical, the Rota-Rack® comes in two versions for different tube sizes. The four colored elements (blue, green, pink, yellow) are connected to each other, but each element can be rotated in the desired direction so that tubes of different sizes can be stored and transported in one rack. Autoclavable.

Rota-Rack® for large tubes

2 x 50 ml

8 x 15 ml

6 x 16 to 20 mm

10 x 10 to 13 mm

Dimensions: 224 x 112 x 112 mm

Additional Information

No.	2-6617
Manufacturer (Brand)	neoLab
EAN	4058072008130
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue yellow green pink
Material	Polypropylene (PP)
Wide	112 mm
Height	112 mm
suitable for	variable applications
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
for vessels with volume up to	15 ml 50 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Combi rack for different tubes and volumes centrifuge tube rack
for vessels with DM up to	20.0 mm 13.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	PCR Reaction Vessels Reaction Vessels Centrifuge Tubes

