




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® SlipLock Flipper® Rack for PCR and reaction tubes, set

Special Price

€44.90 was

€56.50

**plus VAT &
Shipping**

Product Images



Description

The set consists of 4 different colored racks as well as 2 lids that fit over 2 racks each. The racks can be firmly connected to each other by simply sliding them in opposite directions. Can be extended to form a workstation of any length. Stackable.

Each rack provides space for:

- 24 reaction tubes 1.5 or 2.0 ml, max. 11 mm Ø
- 16 reaction tubes 0.5 ml, max. 8 mm Ø
- 60 PCR tubes 0.1 or 0.2 ml
- 5 12-well PCR strips

2 connected racks are ideal for 1 96-well MTP (unskirted, semi-skirted) or 12 8-well PCR strips.

Material: PP, temperature resistant from -80°C to +121°C

Dimensions (L x W x H):

130 x 54 x 35 mm (rack)

133 x 107 x 56 mm (2 racks connected, with lid)

Each set contains 1 rack each in red, blue, green and yellow as well as 2 transparent slip-on lids.

Additional Information

No.	2-2768
Manufacturer (Brand)	neoLab
old neoLab article no.	2-2769
EAN	4058072068967
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue yellow green red
Material	Polypropylene (PP)
Wide	107 mm
Height	56 mm
suitable for	variable applications
Racks features	stackable/collapsible with special chemical resistance with closure/lid available in special sizes for different vessel types/volumes
Number of floors	1
for vessels with volume up to	0.1 ml 0.2 ml 0.5 ml 1.5 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 8.0 mm 11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	PCR reaction vessels test tubes tubes reaction vessels centrifuge tubes

