



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® multi-purpose rack, holding tray + 4 tube racks, neon colored

€79.70
plus VAT &
Shipping

Product Images



Description

Patented rack system consisting of a transparent holding tray and four rack inserts which have different holes on each of the 4 sides:

2 x 30 mm Ø, e.g. for Falcon tubes

6 x 15 mm Ø, e.g. for test tubes 14 mm

9 x 11 mm Ø, e.g. for vials

12 x 10 mm Ø, e.g. for reaction tubes

Material: Polypropylene

Autoclavable up to 121°C

Additional Information

No.	2-2679
Manufacturer (Brand)	neoLab
EAN	4058072068424
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	pink/yellow/green/orange (neon)
Material	Polypropylene (PP)
Wide	10 mm
Height	10 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Centrifuge tube rack test tube rack
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes Reaction vessels Centrifuge tubes
for type vial/vial	Vial. 12 mm Ø

