



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine multifunction rack for tubes with diameter 30 mm, 18 places, yellow

**€14.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

neoLabLine **multifunction rack** – **the rack with many possibilities.**

Suitable for **tubes** with a diameter of 12 to 30 mm or **jars** with a diameter of max. 50 mm.  
Depending on the equipment, a different number of spaces is available.

The positions of the individual compartments are designed so that tubes with screw caps fit perfectly next to each other without interfering. There are matching holes at the bottom of the rack for each tube position, so that all tubes stand securely in all positions.

The height of the compartments can be adjusted and each compartment is alphanumerically coded.

The coding can be continued over 2 connected racks. This connection is made via clamps (accessories).

The **multifunctional rack** can also be used with only one grid instead of 2 (then you need a holding frame as accessory for stabilization)

The multifunction rack is also available as a **basket** (Art. No. 6-2610/6-2611).

With the accessories handle, holding frame and clamps (Art.-No. 6-2612/6-2613) a rack can be converted into a basket.

- Stackable
- Autoclavable

## Additional Information

No.	6-2607
Manufacturer (Brand)	neoLabLine
EAN	4058072156190
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	120 mm
Height	70 mm
Depth	230 mm
Racks features	for different vessel types/volumes alphanumerically coded
Capacity	18
Squares (WxL)	3 x 6
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	30.0 mm
for vessels with volume up to	30 ml
for micro vessel/tube shape	for all shapes
suitable for cooling in	Refrigerator freezer/ freezer

