



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoRack® tube rack interlocking blue/red, PP, 4 x 4 tubes 30 mm

**€10.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The interlocking neoRack® is a practical tube rack that is quick and easy to assemble and disassemble. When disassembled, the neoRack® can be easily stored in a cabinet or drawer in your laboratory. Made of durable polypropylene, this tube rack is temperature resistant from -90°C to +120°C and autoclavable.

The neoRacks are open at the front and are pierced by three rack tiers or square grids that slide into the side and back panels and snap into place. The different models of this neoRack have the same external dimensions 146 x 146 x 70 mm (L x W x H) and differ in the size of the grid. The side and rear panels as well as the bottom shelf are blue for all sizes. The two upper tiers for fixing the tubes are color coded depending on the size.

The easy assembly and disassembly, as well as the stackability make this neoRack® a versatile tube rack for use in your laboratory.

## Additional Information

No.	2-1634
Manufacturer (Brand)	neoRack
EAN	4058072057787
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Wide	146 mm
Height	70 mm
suitable for	variable applications
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
suitable for cooling in	Freezer/ freezer
Tray size	30 x 30 mm
Squares (WxL)	4 x 4
Type frame	Centrifuge tube rack test tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes centrifuge tubes

