



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® vial rack for 50 vials 11.6 mm Ø

€72.80
plus VAT &
Shipping

Product Images



Description

Practical bottle stands made of transparent acrylic, facilitate the transport and handling of sample vials. Safe stand, stackable due to silicone feet. Bottle stands for 8 and 11 mm Ø also suitable for conical vials.

Order No. 7-0649

suitable for: all vials up to 11.6 mm Ø

Threaded vials ND8

Short threaded vials ND9

Threaded vials ND10

Rolled rim vials ND11

Snap ring vials ND11

Additional Information

| | |
|-----------------------------|---|
| No. | 7-0649 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072029142 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | transparent (clear) |
| Material | Polymethylmetacrylate (PMMA) (acrylic glass) |
| Wide | 95 mm |
| Height | 20 mm |
| Racks features | with special chemical resistance suitable for incubators for different vessel types/volumes |
| Number of floors | 1 |
| DM Bore | 12 mm |
| suitable for cooling in | Freezer/ freezer |
| Squares (WxL) | 5 x 10 |
| Type frame | Vial/Vial Rack |
| for vessels with DM up to | 12.0 mm |
| for micro vessel/tube shape | conical. with stand ring round |

