



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® vial rack for 25 vials 24 mm Ø

**€72.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Practical bottle stands made of transparent acrylic**, facilitate the transport and handling of sample vials. Safe stand, stackable due to silicone feet. Bottle stands for 8 and 11 mm Ø also suitable for conical vials.

**Order No. 7-0787**

suitable for:Headspace vials up to 23 mm Ø

Fine thread vials ND18

## Additional Information

---

No.	7-0787
Manufacturer (Brand)	neoLab
EAN	4058072030339
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	160 mm
Height	30 mm
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
DM Bore	24 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	5 x 5
Type frame	Vial/Vial Rack
for vessels with DM up to	24.0 mm
for micro vessel/tube shape	conical. with stand ring round

