



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® BagRack, stackable rack for filter bags

€98.90
plus VAT &
Shipping

Product Images



Description

The **stainless steel rack** is stackable and autoclavable.

Methodical and hygienic storage of bags that have been closed with the ^{BagClip®} is thus guaranteed.

New: Working on opened bags in the rack is possible at all positions.

For up to 12 bags with a width of 19 cm.

Dimensions: 21 x 20 x 36 cm

Additional Information

No.	2-1384
Manufacturer (Brand)	neoLab
EAN	4058072056360
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Wide	205 mm
Height	200 mm
for type bag	BagFilter 400 ml for BagMixer BagPage 400 ml for BagMixer Sterile bags for Bagmixer 400VW

