



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rack for cuvettes 10 mm, 3 x 4 places, model B

€43.60
plus VAT &
Shipping

Product Images



Description

Cuvette rack made of PP for a secure hold of spectrophotometer cuvettes.
The design guarantees a minimum of direct contact.
Only the corners of the cuvettes come into contact with the rack.
Suitable for cuvettes with 10 mm path length.

Order No. 6-2231

Has additionally shaped sides for safe stacking of several racks.
Dimensions: 114 x 111 x 111 mm

Autoclavable up to 121°C

Additional Information

No.	6-2231
Manufacturer (Brand)	neoLab
EAN	4058072025557
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	111 mm
Height	111 mm
Length	114 mm
Depth	114 mm
suitable for	Cuvettes
Racks features	stackable/pluggable with special chemical resistance
suitable for cooling in	Freezer/ freezer
Type frame	Cell rack

