



neoLab® thermometer rack

Product Images



Description

Rack for up to 25 thermometers.
Hole Øtop and middle: 8 mm, bottom: 3.2 mm.
Distance between upper and lower deck: 105 mm

Dimensions: 149 x 213 x 251 mm
Autoclavable up to 121°C

Additional Information

No.	6-2276
Manufacturer (Brand)	neoLab
EAN	4058072025632
Color	Weiß
Material	Polypropylen (PP)
Wide	149 mm
Height	251 mm
Length	213 mm
Depth	213 mm
suitable for	Thermometer
DM Bore	8 mm
Type frame	Thermometergestell

