



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel stand

Special Price

€23.20 was

€34.80

**plus VAT &
Shipping**

Product Images



Description

Stainless steel stand for test tubes. This is a stand in which the tubes automatically align themselves and thus do not touch each other at the "edges" when standing in the stand.

On the one hand, this facilitates the handling of the test tubes and, on the other hand, reduces the risk of contamination from tube to tube.

Autoclavable and disinfectable.

An important aspect in microbiology.

Additional Information

No.	7-3052
Manufacturer (Brand)	neoLab
EAN	4058072077464
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Wide	202 mm
Height	76 mm
Depth	113 mm
suitable for	Test tubes/tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	15.5 mm
suitable for cooling in	Nitrogen tank
Squares (WxL)	5 x 10
Type frame	Test tube rack
for vessels with DM up to	15.5 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

