



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® test tube rack, 31 places, Ø: 16 mm,  
autoclavable, 276x65x110 mm, white, 1 pc**

**Price on  
request**

**plus VAT &  
Shipping**

## Product Images



## Description

---

The white rack for test tubes of different sizes. The places are alphanumerically marked for easy sorting. Made of sturdy polypropylene, and thus autoclavable and chemically resistant to many chemicals.

## Additional Information

---

No.	4-0058
Manufacturer (Brand)	Moonlab
EAN	4058072085766
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	276 mm
Height	110 mm
Depth	65 mm
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	16 mm
Type frame	Test tube rack
for vessels with DM up to	15.5 mm
for vessels with volume up to	16 ml
for micro vessel/tube shape	for all shapes
suitable for cooling in	Refrigerator freezer/ freezer

