



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Moonlab® test tube rack, 20 places, Ø: 20 mm,
autoclavable, 325x74x126 mm, white, 1 pc**

**Price on
request**

**plus VAT &
Shipping**

Product Images



Description

Colorful racks for test tubes of different sizes. The places are alphanumerically marked for easy sorting. Made of sturdy polypropylene, and thus autoclavable and chemically resistant to many chemicals.

Order No. 4-0077: 20 places, white

Additional Information

No.	4-0077
Manufacturer (Brand)	Moonlab
EAN	4058072085919
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	280 mm
Height	126 mm
Depth	74 mm
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	20 mm
for vessels with volume up to	25 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	2 x 10
Type frame	Test tube rack
for vessels with DM up to	19.5 mm
for micro vessel/tube shape	for all shapes

