



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® test tube rack, 12 places, Ø: 32 mm,  
autoclavable, 280x85x126 mm, yellow, 1 pc**

**Price on  
request  
plus VAT &  
Shipping**

## Product Images



## Description

Colorful racks for test tubes of different sizes. The places are alphanumerically marked for easy sorting. Made of sturdy polypropylene, and thus autoclavable and chemically resistant to many chemicals.

**Order No. 4-0082:** 12 places, yellow

## Additional Information

No.	4-0082
Manufacturer (Brand)	Moonlab
EAN	4058072086541
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	280 mm
Height	127 mm
Depth	84 mm
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	32 mm
for vessels with volume up to	50 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	2 x 6
Type frame	Test tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes

