



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Moonlab® test tube rack, 31 places, Ø: 13 mm,
PP, autoclavable, 276x65x110 mm, blue 1 piece**

€9.90
plus VAT &
Shipping

Product Images



Description

Colorful racks for test tubes of different sizes. The places are alphanumerically marked for easy sorting. Made of sturdy polypropylene, and therefore autoclavable and chemically resistant to many chemicals.

Order No. 4-0090: 31 places, blue

Additional Information

No.	4-0090
Manufacturer (Brand)	Moonlab
EAN	4058072087357
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Wide	276 mm
Height	110 mm
Depth	65 mm
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	13 mm
for vessels with volume up to	10 ml
suitable for cooling in	Refrigerator freezer/ freezer
Type frame	Test tube rack
for vessels with DM up to	12.5 mm
for micro vessel/tube shape	for all shapes

